

Double Loop High Back Sling



Operating Instructions

To Reorder

Wy'East Medical Corporation
P.O. Box 1625
16700 S.E. 120th Ave.
Clackamas, OR 97015, USA

(503) 657-3101

(800) 255-3126

FAX: (503) 657-6901

www.wyeastmed.com
contact@wyeastmed.com

For more detailed instructions on how to properly operate and use a sling with your lift, please refer to the operating manual provided with your lift.

WY'EAST Medicals policy is one of continuous development and improvement, and therefore, reserves the right to change these specifications without notice.



Read carefully the safety instructions and warnings before attempting to use the Wy'East Medical DOUBLE LOOP HIGH BACK SLING.

SAFETY INSTRUCTIONS AND WARNINGS

These instructions are to be used in conjunction with those issued with your vertical lift system.

□ Inspect the sling carefully before each use (see figure 1). If fraying, seam separation, or other damage is visible or suspected, dispose of it immediately and replace it with a new DOUBLE LOOP HIGH BACK SLING.

□ The Wy'East Medical DOUBLE LOOP HIGH BACK SLING is rated up to a maximum weight capacity of 600 lbs (273 kg) & up to 800 lbs (362 kg) for the X-Large Sling.

□ Never exceed the weight limitation of your lift or sling, as serious harm to the patient or caregiver may result.

STANDARD REUSABLE SLINGS

LAUNDRY INSTRUCTIONS

- Machine wash in warm water and mild detergent,
- Warm tumble-dry.

Wy'East Medical DOUBLE LOOP HIGH BACK SLINGS have been designed to be used in a manner that is consistent with that of your vertical lift system.

Wy'East Medical DOUBLE LOOP HIGH BACK SLINGS are to be used only with vertical lift systems that use a sling lifting bar with a 2 point attachment as shown in figure 3.

Before using the DOUBLE LOOP HIGH BACK SLING, familiarize yourself with the proper method of sling operation found in the operator's manual provided with your lift system.

The patient to be transferred must be assessed by a qualified professional before attempting to transfer.


Always inspect the sling before every use for any wear or damage and replace as needed (see figure 1).

Place the sling under the patient by using a rolling technique, if the patient is on the floor or in a bed. If the patient is sitting up in a chair, slide the sling down the patient's back from behind the patient's head.

Bring the leg bands under and up between the thighs.

Route the leg bands through each other as shown in figure 2 and connect to opposite ends of the sling lifting bar.

Connect the loops at head end of sling to the same side of sling lifting bar as shown in figure 3.



**FIGURE 1
SLING INSPECTION POINTS**



**FIGURE 2
CROSSED LEG BANDS**



**FIGURE 3
LOOP CONNECTION**