



The TC-300, with its 400-lb. transport capacity, converts from a stretcher to a chair in seconds.



The unique design of the TC-300 offers the ultimate in patient transfer flexibility. Adherence to the highest manufacturing standards gives you a product you'll depend on for years.



The TC-300 with optional Stirrup/Leg Drop and two-piece mattress pad allows the clinician to perform pelvic exams.


TC-300 Specifications

The adjustable positioning of the TC-300 offers the ultimate in transport flexibility.

- HIGH POSITION:37"
- LOW POSITION:23"
- WIDTH OF PAD:23"
- MAXIMUM WIDTH:27"
- MAX. PATIENT WEIGHT:400 LBS.

The TC-300

The all-purpose treatment chair that converts into a stretcher



The TC-300, an all-purpose treatment chair that easily converts into a stretcher, adapts instantly to your departmental needs for patient examination, treatment, and transport. Designed for maximum patient comfort and absolute convenience for your clinical staff, the TC-300 will bring ease and security to every procedure.

- **COMPLETELY ADJUSTABLE.** The TC-300 quickly adapts to accommodate procedure or patient need from flat to upright positions. Height adjustment ranges from 23" to 37" and is controlled by a convenient foot pedal. And the Trendelenburg activates quickly – in one easy step – causing no delays in emergency situations.
- **EASY TO USE.** The broad wheelbase ensures stability, but still leaves plenty of room to maneuver through small rooms and narrow or crowded halls. Controls are placed on both sides for easy accessibility during procedures or transport.
- **OPTIONAL ACCESSORIES.** Accessories include stirrup/leg drop and two-piece mattress pad, tray, upper side rails, padded cradle arm supports, head rest, and comfort pad mattress.

Learn more about the TC-300 today. Call 1-800-255-3126 or visit www.wyeastmed.com.

