





TOTALIFT-II & TC-300 (All Models) Recycling Instructions



TOTALIFT-II & TC-300

Options


Recycling Instructions



5-00990-00
Upper Side-Rail



Metals



51124
Leg-Drop/Stirrup



Metals



5-01023-00
Patient Tray



Metals
Plastics



50841
IV Pole



Metals



5-00932-00
Safety Strap



Plastics



5-01027-00
Arm Cradle



Metals

Mattress (All Models) – Combustible waste

Headrest (All Models) – Combustible waste